


Flex™ Product Portfolio



Enhanced Beating Heart and Surgical Stabilization

The new Terumo Flex beating heart and surgical stabilization product portfolio delivers versatile performance and exceptional control during cardiac surgery.

Terumo Flex Product Portfolio

Beating Heart Products

Herc Flex™ Stabilizer Arm

The Herc Flex stabilizer arm provides a multifunctional platform for both on- and off-pump cardiac surgery support. Its lower-profile, smaller diameter and lighter weight design allows precise positioning, stabilization and retraction for exceptional visualization in the surgical field. The reengineered handle improves control when tightening and maintains maximum stability. The Herc Flex stabilizer arm supports off-pump disposables, including Atlas™ Positioner Attachments, Titan and Titan Flex Stabilizers, and Terumo Assistant™ Attachment disposable products.



Herc Flex™ Stabilizer Arm

Titan™ Flex Stabilizer

The Titan Flex stabilizer is redesigned to deliver greater versatility and articulation during off-pump cardiac surgery. It features an enhanced malleable framework that contours to a variety of heart anatomy and maintains excellent stability for exposure and anastomosis.



Titan™ Flex Stabilizer

Ordering Information

Catalog #	Description	Units/Case
Herc Flex Stabilizer Arm		
T401162	Herc Flex 360 Stabilizer Arm, universal; Top mount, low profile, 360-degree rotation, reusable stainless steel	1
T401163	Herc Flex Stabilizer Arm, universal; Top mount, low profile, reusable stainless steel	1
Titan Flex Stabilizer		
T401232U	Titan Flex 360 Stabilizer; Single use, suction, 360-degree articulation maximizes vessel and stabilization, enhanced malleable footprint	1
T401232	Titan Flex Stabilizer; Single use, suction, enhanced malleable footprint	1

Surgical Stabilization Products

Assistant Flex Attachment with StableSoft™ Technology

The Assistant Flex allows more custom configurations compared to existing Assistant products, with enhanced malleable zones at the base of the device for precise contouring. The textured, silicone-coated surface provides exceptional grip and stability and is ideal for surgical cases requiring an extra hand.



Assistant™ Flex Attachment

Ordering Information

Catalog #	Description	Units/Case
Assistant Flex Attachment with StableSoft		
T401213S	Assistant Flex Attachment with StableSoft Technology, 3L; Single use, fully malleable, three fingers, long	5
T401214S	Assistant Flex Attachment with StableSoft Technology, 3S; Single use, fully malleable, three fingers, short	5



Terumo Cardiovascular Group

6200 Jackson Road, Ann Arbor, MI 48103-9300 USA
Tel: +1.734.663.4145 Fax: +1.734.663.7981
Email: cardiovascular@terumomedical.com

For a complete list of offices, please consult our website.

www.terumo-cvgroup.com/offices

Terumo Corporation

2-44-1 Hatagaya, Shibuya-Ku, Tokyo 151-0072 Japan
Tel: +81.3.3374.8111 Fax: +81.3.3374.8399

Terumo Europe NV

Authorized EC Representative
Interleuvenlaan 40, 3001 Leuven, Belgium
Tel: +32.16.38.12.11 Fax: +32.16.40.02.49

Terumo Asia Holdings Pte. Ltd.

300 Beach Road, #33-03, The Concourse Singapore 199555
Tel: +65.6.295.1792 Fax: +65.6.294.2329

Terumo Latin America Corporation

8750 NW 36th Street, Suite 600, Miami, FL 33178 USA
Tel: +1.305.477.4822 Fax: +1.305.477.4872